# BSNScreen Card fFN Test

Reagenti per 10 test

Manuale d'uso

**REF**

TR1002-10

**IVD**

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**CE**

**B.S.N. BIOLOGICAL SALES NETWORK S.R.L.**  
26012 Castelleone (CR) Italia - Via Coelli, 18  
Tel. +39 0374 351005 - Fax +39 0374 57965 - e-mail: [info@bsn-srl.it](mailto:info@bsn-srl.it) –  
[bsn@postecert.it](mailto:bsn@postecert.it)  
VAT: 11317290150 - R.E.A. di Cremona n. 143395

**TR1002-10**

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## 9 BIBLIOGRAFIA

1. Friedman ML, McElin TW. *Diagnosis of ruptured fetal membranes. Am J Obstet Gynecol* 1969;104:544-550
2. Steinman G, Kleiner GJ, Greston WM. *Spontaneous rupture of membranes. NY St. J Med* 1979;1849-1851
3. Koninckx PR, Trappeniers H, van Assche FA. *Prolactin concentration in vaginal fluid: a new method for diagnosing ruptured membranes. Br J Obstet Gynaecol* 1981;88:607-610
4. Rochelson BL, Richardson DA, Macri JN. *Rapid assay: possible application in the diagnosis of premature rupture of the membranes. Obstet Gynecol* 1983;62:414-418
5. Bell SC, Hales MW, Patel SR, Kirwan PH, Drife JO, Milford-Ward A. *Amniotic fluid concentrations of secreted pregnancy-associated endometrial alpha1 and alpha2 globulins (alpha1 and alpha2 PEG). Br J Obstet Gynaecol* 1986;93:909-915
6. Suikkari AM, Rutanen EM, Seppälä M. *Circulating levels of immunoreactive insulin-like growth factor-binding protein in non-pregnant women. Hum Reprod* 1987;2:297-300
7. Bell SC, Keyte JW; *N-terminal amino acid sequence of human pregnancy-associated endometrial alpha1-globulin, an endometrial insulin-like growth factor (IGF) binding protein – Evidence for two small molecular weight IGF binding proteins. The Endoc Soc* 1988;1202-1204
8. Rutanen EM, Pekonen F, Kärkkäinen T. *Measurement of insulin-like growth factor binding protein-1 in cervical/vaginal secretions: comparison with the ROM-check membrane immunoassay in the diagnosis of ruptured fetal membranes. Clin Chim Acta* 1993;214:73-81
9. Medina TM, Hill DA. *Preterm premature rupture of membranes: diagnosis and management. Am Fam Physic* 2006;73: 659-664
10. Lockwood CJ, Wein R, Chien D, Ghidini A, Alvarez M, Berkowitz RL. *Fetal membrane rupture is associated with the presence of insulin-like growth factor-binding protein-1 in vaginal secretions. Am J Obstet Gynecol* 1994;171:146-150
11. Guibourdenche J, Luton D, André E, Noël M, Porquet D. *Détection rapide de l'Insulin-like growth factor-binding protein-1 et de la fibronectine fœtale dans les sécrétions cervico-vaginales pour le diagnostic de la rupture prématurée des membranes. Ann Clin Biochem* 1999 ;36 :388-390
12. Sucak A, Moroy P, Cakmak P, Mungan T, Mollamahmutoglu L, Danisman N. *Insulin-Like Growth Factor Binding Protein-1; a rapid detection of amniotic fluid leakage after amniocentesis. Turk J Med Sci* 2005;35:157-161
13. Talabani H, Recoules A, Grenier S, Luton D, Paper T, Porquet D, Oury JF, Guibourdenche J. *L'alpha1-microglobuline placentaire dans les sécrétions vaginales : un nouveau marqueur de la rupture prématurée des membranes. Act Pharm Biol Clin* 2005 :12 :233-240
14. Akercan F, Cirpan T, Kazandi M, Terek MC, Mgoyi I, Ozkinay E. *The value of insulin-like growth factor binding protein-1 in the cervical-vaginal secretion detected by immunochromatographic dipstick test in the prediction of delivery in women with clinically unconfirmed preterm premature rupture of membranes. Eur J Obstet Gynecol Rep Biol* 2005;121:159-163

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