

Order No.	Description	Quantity
EUM08100	FloMass® Drug of Abuse in Whole Blood (Amphetamine, Methamphetamine, MDA, MDE, MDMA, MBDB, BEG, Cocaine, EME, Cocaethylene, Morphine, Codeine, Dihydrocodeine, 6-MAM, THC, 11-OH-THC, THC-COOH, Methadone, EDDP, Buprenorphine, Norbuprenorphine, Ketamina) LC-MS/MS detection	100 assays
	Contents	
	Mobile Phase A	1 x EUM02011
	Mobile Phase B	1 x EUM02012
	Solution 1 Buffer	1 x EUM08021
	Solution 2 Precipitant	1 x EUM08022
	Internal Standard Mix	1 x EUM08031
	Separately available components	
EUM02011	Mobile Phase A	600 mL
EUM02012	Mobile Phase B	500 mL
EUM08021	Solution 1 Buffer	12 mL
EUM08022	Solution 2 Precipitant	65 mL
EUM08031	Internal Standard Mix	2x1 mL
	Accessory	
EUM00C02	Analytical Column with test chromatogram	1 pc
EUM00A04	Precolumn	4 pcs
EUM00A05	Holder (incl. 1 precolumn)	1 pc
EUM02032	Internal Standard Mix Conc.	2x0.75 mL
EUM02071	FloTuning Mix A + Mix B	2x2 mL
12WH0720A	Control Set for Drug of Abuse in Whole Blood, lyoph.	10x2.5 mL
12WH0618B	Control Set for Drug of Abuse in Whole Blood, lyoph.	10x2.5 mL
EUM08055	Control Set for Drug of Abuse in Whole Blood, lyoph. (2 levels)	3x2x0.6 mL
EUM08042	Calibrator Set for Drug of Abuse in Whole Blood, lyoph. (6 levels)	3x6x0.6 mL