

Order No.	Description	Quantity
<b>EUM14200</b>	<b>FloMass® Cortisol/Cortisone in Saliva</b> LC-MS/MS detection	200 assays
	<b>Contents</b>	
	Mobile Phase A	1 x EUM06011
	Mobile Phase B	1 x EUM06012
	Solution 1 Dilution	1 x EUM14021
	Solution 2 Extraction	1 x EUM14022
	Solution 3 Reconstitution	1 x EUM14023
	Internal Standard	1 x EUM14031
	<b>Separately available components</b>	
EUM06011	Mobile Phase A	900 mL
EUM06012	Mobile Phase B	900 mL
EUM14021	Solution 1 Dilution	60 mL
EUM14022	Solution 2 Extraction	250 mL
EUM14023	Solution 3 Reconstitution	10 mL
EUM14031	Internal Standard	4.5 mL
	<b>Accessory</b>	
EUM00C14	Analytical Column with test chromatogram	1 pc
EUM00A06	Precolumns	4 pcs
EUM00A07	Holder (incl. 1 precolumn)	1 pc
EUM06071	FloTuning	1x2 mL
EUM14051	Control Set for Cortisol/Cortisone in Saliva (3 levels)	2x3x2 mL
EUM14041	Calibrator Set for Cortisol/Cortisone in Saliva, lyoph. (6 levels)	2x6x1 mL