

Order No.	Description	Quantity
EUM17100	FloMass® Catecholamines, Free Metanefrines and Serotonin in Urine LC-MS/MS detection	100 assays
	Contents	
	Mobile Phase A	1 x EUM17011
	Mobile Phase B	1 x EUM17012
	Complexing Reagent	1 x EUM17021
	pH Correcting Reagent	1 x EUM17022
	SPE Activating Reagent	1 x EUM17023
	SPE Conditioning Reagent	1 x EUM17024
	SPE Washing Reagent	1 x EUM17025
	SPE Eluting Reagent	1 x EUM17026
	Internal Standard Mix	1 x EUM17031
	SPE Columns	1 x EUM17061
	Separately available components	
EUM17011	Mobile Phase A	500 mL
EUM17012	Mobile Phase B	500 mL
EUM17021	Complexing Reagent	50 mL
EUM17022	pH Correcting Reagent	10 mL
EUM17023	SPE Activating Reagent	110 mL
EUM17024	SPE Conditioning Reagent	110 mL
EUM17025	SPE Washing Reagent	110 mL
EUM17026	SPE Eluting Reagent	65 mL
EUM17031	Internal Standard Mix	2x1 mL
EUM17061	SPE Columns	100 pcs
	Accessory	
EUM00C17	Analytical Column with test chromatogram	1 pc
EUM00A14	Precolumns	5 pcs
EUM00A15	Holder (incl. 1 precolumn)	1 pc
EUM17071	FloTuning Mix A + Mix B	2x2 mL
EUM17051	Control Set for Biogenic Amines in Urine, lyoph. (Adrenaline, Noradrenaline, Dopamine, Metanephrine, Normetanephrine, 3-Methoxytyramine, VMA, HVA, 5-HIAA and Creatinine) (2 levels)	2x2x2 mL
EUM17041	Calibrator Set for Catecholamines, Metanefrines and Serotonin in Urine in Urine, lyoph. (7 levels)	2x7x1 mL

Vers.: 01/07/2022