

Order No.	Description	Quantity
<b>EUH01100</b>	<b>FloChrom® Catecholamines in Urine</b> Fluorimetric detection	100 assays
	<b>Contents</b>	
	Mobile Phase	1 x EUH01011
	Reagent 1	1 x EUH01021
	Reagent 2	1 x EUH01022
	Protective Reductant	1 x EUH01023
	Internal Standard	1 x EUH01031
	Calibration	1 x EUH01041
	SPE Column	1 x EUH01061
	<b>Separately available components</b>	
EUH01011	Mobile Phase	2000 mL
EUH01021	Reagent 1	200 mL
EUH01022	Reagent 2	60 mL
EUH01023	Protective Reductant	7 g
EUH01031	Internal Standard	10 mL
EUH01041	Calibration	10x2.5 mL
EUH01061	SPE Column	100 pcs
	<b>Accessory</b>	
EUH01090	Analytical Column with test chromatogram	1 pc
EUH01070	Precolumns	5 pcs
EUH01080	Holder (incl. 1 precolumn)	1 pc
	<b>ClinChek® Controls</b>	
8820	Urine Control, lyophil. (Level I)	10x8 mL
8821	Urine Control, lyophil. (Level II)	10x8 mL
8822	Urine Control, lyophil. (Level I, II)	2x5x8 mL